

The Duraline Recirculating Lab Faucet comes with the Recirculating Module and Tube temporarily assembled.

**STEP 1.** Remove (unscrew) the Recirculating Module and Recirculating Tube from the through deck nipple.

**STEP 2.** Mount the faucet in the desired location through a 1" diameter hole. Tighten to deck with the Ring Nut supplied.

**STEP 3.** Insert the Recirculating Tube into and through the Through Deck Nipple until the 3/8" male NPT thread mates with the 3/8" female NPT thread in the inlet connection of the Through Deck Nipple. Tighten by hand plus 1/4 wrench turn.


**STEP 4.** Connect a 3/8" supply line tube to the 3/8" inlet compression tube fitting.

**STEP 5.** Connect a 1/4" return line tube to the 1/4" outlet compression tube fitting.

**STEP 6.** Pressure test all connections.

Note: A high density PTFE tape is the recommended thread sealant.

SECTION A-A  
SCALE 1 : 4

Designer			TITLE: Duraline Recirculating Lab Faucet, Deck Mount, With Recirculating Module		
Review					
Engineer					
			DWG. NO.	DL-DDR-2 (Natural Polypropylene) DL-DDR-3 (Natural PVDF)	REV
			UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES	2950 Airway Avenue, Costa Mesa, CA 92626	SCALE: N/A